



# SLPI Antibody Pair

<b>Product Code</b>	CSB-EAP07154
<b>Species Reactivity</b>	Human,Bovine,Horse,Pig,Chicken,Goat
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	ALP,BLPI,HUSI-1,Mucus proteinase inhibitor,MPI,Protease inhibitor WAP4,Secretory leukocyte protease inhibitor,Seminal proteinase inhibitor,WAP four-disulfide core domain protein 4,SLPI,WAP4, WFDC4
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Gene Names</b>	SLPI
<b>Protein Names</b>	Antileukoproteinase
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP07154C Detection: CSB-EAP07154D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4