



PAX2 Antibody Pair

Product Code	CSB-EAP11564
Species Reactivity	Human,Horse,Pig,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	PAX2
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	PAX2
Protein Names	Paired box protein Pax-2
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP11564C Detection: CSB-EAP11564D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4