



# VAMP7 Antibody Pair

<b>Product Code</b>	CSB-EAP025785
<b>Species Reactivity</b>	Human,Bovine,Horse,Pig,Goat
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	Synaptobrevin-like protein 1,Tetanus-insensitive VAMP
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Gene Names</b>	VAMP7
<b>Protein Names</b>	Vesicle-associated membrane protein 7
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP025785C Detection: CSB-EAP025785D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4