



VAMP7 Antibody Pair

Product Code	CSB-EAP025785
Species Reactivity	Human,Bovine,Horse,Pig,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Synaptobrevin-like protein 1,Tetanus-insensitive VAMP
Product Type	Antibody Pairs
Species	Homo sapiens (Human)
Gene Names	VAMP7
Protein Names	Vesicle-associated membrane protein 7
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP025785C Detection: CSB-EAP025785D
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4