



PTK2 Antibody Pair

Product Code	CSB-EAP14445
Species Reactivity	Guinea pig,Bovine,Horse,Pig,Dog,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Focal adhesion kinase-related nonkinase,p125FAK,pp125FAK
Product Type	Antibody Pairs
Immunogen Species	Gallus gallus (Chicken)
Gene Names	PTK2
Protein Names	Focal adhesion kinase 1
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP14445C Detection: CSB-EAP14445D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4