



KCNA1 Antibody Pair

Product Code	CSB-EAP012005
Species Reactivity	Bovine,Horse,Dog,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Voltage-gated K(+) channel HuK1,Voltage-gated potassium channel HBK1,Voltage-gated potassium channel subunit Kv1.1,KCNA1
Product Type	Antibody Pairs
Species	Homo sapiens (Human)
Gene Names	KCNA1
Protein Names	Potassium voltage-gated channel subfamily A member 1
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP012005C Detection: CSB-EAP012005D
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4