



Gstp1 Antibody Pair

| | |
|----------------------------|---|
| Product Code | CSB-EAP009989 |
| Species Reactivity | Rat,Bovine,Horse,Pig |
| Tested Applications | S-ELISA |
| Form | Liquid |
| Alias | Chain 7,GST 7-7,GST class-pi,Gstp1 |
| Product Type | Antibody Pairs |
| Immunogen Species | Rattus norvegicus (Rat) |
| Target Names | Gstp1 |
| Protein Names | Glutathione S-transferase P |
| Notes | We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.125ug/ml.Optimal dilutions should be determined experimentally by the researcher. |
| Host | Capture: Rabbit Detection: Rabbit |
| Components | Capture: CSB-EAP009989C Detection: CSB-EAP009989D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol. |
| Storage-Buffer | Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |