



XRCC6 Antibody Pair

Product Code	CSB-EAP01617
Species Reactivity	Bovine,Horse
Tested Applications	S-ELISA
Form	Liquid
Alias	5'-deoxyribose-5-phosphate lyase Ku70,70 kDa subunit of Ku antigen,ATP-dependent DNA helicase 2 subunit 1,ATP-dependent DNA helicase II 70 kDa subunit ,CTC box-binding factor 75 kDa subunit
Product Type	Antibody Pairs
Species	Homo sapiens (Human)
Gene Names	XRCC6
Protein Names	X-ray repair cross-complementing protein 6 ptein
Notes	We recommend using the capture antibody at a concentration of 0.1ug/ml and the detection antibody at a concentration of 0.1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01617C Detection: CSB-EAP01617D
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4