



CREB1 Antibody Pair

Product Code	CSB-EAP005947
Species Reactivity	Human,Guinea pig,Bovine,Horse,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	CREB1
Product Type	Antibody Pairs
Species	Homo sapiens (Human)
Gene Names	CREB1
Protein Names	Cyclic AMP-responsive element-binding protein 1
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP005947C Detection: CSB-EAP005947D
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4