



L1CAM Antibody Pair

Product Code	CSB-EAP07439
Species Reactivity	Human,Bovine,Horse,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	INCAM-100,CD_antigen=CD106 VCAM1,L1CAM
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	L1CAM
Protein Names	Vascular cell adhesion protein 1
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07439C Detection: CSB-EAP07439D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4