



ESR1 Antibody Pair

Product Code	CSB-EAP11399
Species Reactivity	Human,Horse,Pig,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	ER-alpha,Estradiol receptor,Nuclear receptor subfamily 3 group A member 1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	ESR1
Protein Names	Estrogen receptor
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP11399C Detection: CSB-EAP11399D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.