



IL1B Antibody Pair

Product Code	CSB-EAP011614
Species Reactivity	Guinea pig,Chicken,Pig,Dog,Bovine,Sheep,Goat,Horse
Tested Applications	S-ELISA
Form	Liquid
Alias	IL1B
Product Type	Antibody Pairs
Immunogen Species	Ovis aries (Sheep)
Gene Names	IL1B
Protein Names	Interleukin-1 beta
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP011614C Detection: CSB-EAP011614D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4