



CST3 Antibody Pair

Product Code	CSB-EAP06197
Species Reactivity	Human,Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	Cystatin-3,Gamma-trace,Neuroendocrine basic polypeptide,Post-gamma-globulin,CST3
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	CST3
Protein Names	Cystatin-C
Notes	We recommend using the capture antibody at a concentration of 0.35ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Mouse Detection: Goat
Components	Capture: CSB-EAP06197C Detection: CSB-EAP06197D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4