



Timp1 Antibody Pair

Product Code	CSB-EAP167191
Species Reactivity	Rabbit,Rat
Tested Applications	S-ELISA
Form	Liquid
Alias	Metalloproteinase inhibitor 1,Tissue inhibitor of metalloproteinases 1,TIMP-1,Timp1
Product Type	Antibody Pairs
Immunogen Species	Rattus norvegicus (Rat)
Target Names	Timp1
Protein Names	Metalloproteinase inhibitor 1
Notes	We recommend using the capture antibody at a concentration of 4ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Mouse Detection: Mouse
Components	Capture: CSB-EAP167191C Detection: CSB-EAP167191D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.