



RBP4 Antibody Pair

Product Code	CSB-EAP070501
Species Reactivity	Human,Horse,Dog,Chicken,Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	Plasma retinol-binding protein, PRBP, RBP, RBP4, PRO2222
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RBP4
Protein Names	Retinol-binding protein 4
Notes	We recommend using the capture antibody at a concentration of 0.6ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Mouse Detection: Mouse
Components	Capture: CSB-EAP070501C Detection: CSB-EAP070501D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.