



Vegfc Antibody Pair

Product Code	CSB-EAP07725
Species Reactivity	Human,Mouse,Bovine,Horse,Pig,Dog,Duck
Tested Applications	S-ELISA
Form	Liquid
Alias	Vascular endothelial growth factor C,VEGF-C,Flt4 ligand,Flt4-L,Vascular endothelial growth factor-related protein,VRP,Vegfc
Product Type	Antibody Pairs
Species	Mus musculus (Mouse)
Gene Names	Vegfc
Protein Names	Vascular endothelial growth factor C
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07725C Detection: CSB-EAP07725D
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4