



Vegfc Antibody Pair

| | |
|----------------------------|--|
| Product Code | CSB-EAP07725 |
| Species Reactivity | Human,Mouse,Bovine,Horse,Pig,Dog,Duck |
| Tested Applications | S-ELISA |
| Form | Liquid |
| Alias | Vascular endothelial growth factor C,VEGF-C,Flt4 ligand,Flt4-L,Vascular endothelial growth factor-related protein,VRP,Vegfc |
| Product Type | Antibody Pairs |
| Immunogen Species | Mus musculus (Mouse) |
| Gene Names | Vegfc |
| Protein Names | Vascular endothelial growth factor C |
| Notes | We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher. |
| Host | Capture: Rabbit Detection: Rabbit |
| Components | Capture: CSB-EAP07725C Detection: CSB-EAP07725D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol. |
| Storage-Buffer | Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4 |