



# Recombinant *Saccharomyces cerevisiae* Phosphatidylinositol N- acetylglucosaminyltransferase ERI1 subunit (ERI1)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-CF351787SVG  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P62651   |
| <b>Product Type</b>      | Transmembrane Protein  |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)   |
| <b>Sequence</b>          | MRPRDQGFLVLGFTYSVLLISLATFYWLRNNDNFLHYWCVLLLCPATLWLWALI<br>AWCDSE MFASSKDE  |
| <b>Source</b>            | in vitro E.coli expression system  |
| <b>Target Names</b>      | ERI1   |
| <b>Protein Names</b>     | Recommended name: Phosphatidylinositol N-acetylglucosaminyltransferase<br>ERI1 subunit Alternative name(s): Endoplasmic reticulum-associated Ras<br>inhibitor protein 1  |
| <b>Expression Region</b> | 1-68   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at<br>4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,<br>storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life<br>of lyophilized form is 12 months at -20°C/-80°C. |