



Recombinant Coffea arabica Photosystem II reaction center protein H (psbH)

Product Code	CSB-CF366319DVA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0A364
Product Type	Transmembrane Protein
Immunogen Species	Coffea arabica (Arabian coffee)
Sequence	ATQTVDNSSRSGPRRTTVGSLLKPLNSEYGKVVAPGWGTTPLMGVAMALFAVF LSIILEIY NSSVLLDGISMN
Source	in vitro E.coli expression system
Target Names	psbH
Protein Names	Recommended name: Photosystem II reaction center protein H Short name= PSII-H Alternative name(s): Photosystem II 10 kDa phosphoprotein
Expression Region	2-73
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.