



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein SCRG_04509 (SCRG_04509)

Product Code	CSB-CF459238SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LS52
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	MSKHKHEWTESVANS GPASILSYCASSILMTVTNKFVVNLDNFMNMFVMLFVQ SLVCTVT LCILRIVGVANF
Source	in vitro E.coli expression system
Target Names	SCRG_04509
Protein Names	Recommended name: Uncharacterized protein SCRG_04509
Expression Region	1-72
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.