



Recombinant Pongo pygmaeus Cytochrome c oxidase subunit 6C (COX6C)

Product Code	CSB-CF801780EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7YRK3
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	APEVLPKPRMRGLLARRLRNHMAVAFVLSLGVAALYKFRVADRRKKAYADFY RNYDVMKD FEEMRKAGIFQSVK
Source	in vitro E.coli expression system
Target Names	COX6C
Protein Names	Recommended name: Cytochrome c oxidase subunit 6C Alternative name(s): Cytochrome c oxidase polypeptide VIc
Expression Region	2-75
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.