



Recombinant *Oryza sativa* subsp. *indica* ATP synthase subunit 9, mitochondrial (ATP9)

Product Code	CSB-CF313945OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C518
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MLEGAKLIGAGAATIALAGAAVIGIGNVFSSLIHSVARNPSLAKQLFGYAILGFAL TEAIA LFALMMAFLILFVF
Source	in vitro <i>E.coli</i> expression system
Target Names	ATP9
Protein Names	Recommended name: ATP synthase subunit 9, mitochondrial Alternative name(s): Lipid-binding protein
Expression Region	1-74
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.