



Recombinant Lama guanicoe pacos Cytochrome c oxidase subunit 6C (COX6C)

Product Code	CSB-CF611310LBU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0Q4Z0
Product Type	Transmembrane Protein
Immunogen Species	Vicugna pacos (Alpaca) (Lama pacos)
Sequence	SSGALLPKPQMRGLLAKRLRVHIVGAFVVALGVAAAYKFGVAEPRKKAYADFY RNYDSMK DFEEMRQAGVFQSAK
Source	in vitro E.coli expression system
Target Names	COX6C
Protein Names	Recommended name: Cytochrome c oxidase subunit 6C Alternative name(s): Cytochrome c oxidase polypeptide VIc
Expression Region	2-76
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.