



Recombinant Pongo abelii ATP synthase lipid-binding protein, mitochondrial (ATP5G2)

Product Code	CSB-CF716141PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5RAP9
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	DIDTAAKFIGAGAATVGVAGSGAGIGTVFGSLIIGYARNPSLKQQLFSYAILGFA LSEAM GLFCLMVAFLILFAM
Source	in vitro E.coli expression system
Target Names	ATP5G2
Protein Names	Recommended name: ATP synthase lipid-binding protein, mitochondrial Alternative name(s): ATP synthase proteolipid P2 ATPase protein 9 ATPase subunit c
Expression Region	67-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.