



# Recombinant Pongo abelii ATP synthase lipid-binding protein, mitochondrial (ATP5G3)

<b>Product Code</b>	CSB-CF716241PYX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5RFL2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Sequence</b>	DIDTAAKFIGAGAATVGVAGSGAGIGTVFGSLIIGYARNPSLKQQLFSYAILGFA LSEAM GLFCLMVAFLILFAM
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ATP5G3
<b>Protein Names</b>	Recommended name: ATP synthase lipid-binding protein, mitochondrial Alternative name(s): ATP synthase proteolipid P3 ATPase protein 9 ATPase subunit c
<b>Expression Region</b>	68-142
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.