



Recombinant Pig ATP synthase lipid-binding protein, mitochondrial (ATP5G1)

Product Code	CSB-CF002359PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1XQS5
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	DIDTAAKFIGAGAATVGVAGSGAGIGTVFGSLIIGYARNPSLKQQLFSYAILGFA LFEAM GLFCLMVAFLILFAM
Source	in vitro E.coli expression system
Target Names	ATP5G1
Protein Names	Recommended name: ATP synthase lipid-binding protein, mitochondrial Alternative name(s): ATP synthase proteolipid P1 ATPase protein 9 ATPase subunit c
Expression Region	62-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.