



Recombinant Duck hepatitis B virus Large envelope protein (S)

Product Code	CSB-CF365861DHJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03145
Product Type	Transmembrane Protein
Immunogen Species	Duck hepatitis B virus (strain United States/DHBV-16) (DHBV)
Sequence	MSGTFGGILAGLIGLLVSFFLLIKILEILRRLDWWWISLSSPKGKMQCAFQDTGA QISPH YVGSCPWGCPGFLWTYL
Source	in vitro E.coli expression system
Target Names	S
Protein Names	Recommended name: Large envelope protein Alternative name(s): L glycoprotein L-HBsAg Short name= LHB Large S protein Large surface protein Major surface antigen Cleaved into the following chain: 1. Truncated S protein
Expression Region	162-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.