



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YMR158W-B (YMR158W-B)

Product Code	CSB-CF718356SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6B0X2
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	WECSFFRSESFCKTLFSMVPLINSSFNLSVFLFFNAITSFNLRISCSLLFSSFC RIASV FNKASSWLTMPPMAPLLS
Source	in vitro E.coli expression system
Target Names	YMR158W-B
Protein Names	Recommended name: Putative uncharacterized protein YMR158W-B
Expression Region	28-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.