



Recombinant UPF0197 transmembrane protein Y57E12AM.1 (Y57E12AM.1)

Product Code	CSB-CF839264CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q965T1
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MDISKMNRYTAPVNFASLPLLTFTLCGVGLLLLATFTMIQVTSTKYNRNLLKELF IAATS SVFLGFGSVFLLLWVGIVV
Source	in vitro E.coli expression system
Target Names	Y57E12AM.1
Protein Names	Recommended name: UPF0197 transmembrane protein Y57E12AM.1
Expression Region	1-79
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.