



Recombinant Pongo abelii V-type proton ATPase subunit e 1 (ATP6V0E1)

Product Code	CSB-CF716145PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5RAV0
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	AYHGLTVPLIVMSVFWGFGFLVPWFIPKGPNGVITMLVTCSVCCYLFWLI LAQLN PLFGPQLKNETIWYLKYHWP
Source	in vitro E.coli expression system
Target Names	ATP6V0E1
Protein Names	Recommended name: V-type proton ATPase subunit e 1 Short name= V-ATPase subunit e 1 Alternative name(s): Vacuolar proton pump subunit e 1
Expression Region	2-81
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.