



Recombinant Rat V-type proton ATPase subunit e 1 (Atp6v0e1)

Product Code	CSB-CF755280RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q794C0
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	AYHGLTVPLIVMSVFWGFGVLLVPWFIPKGPNGVITMLVTCSVCCYLFWLIAI LAQLN PLFGPQLKNETIWYLKYHWP
Source	in vitro E.coli expression system
Target Names	Atp6v0e1
Protein Names	Recommended name: V-type proton ATPase subunit e 1 Short name= V-ATPase subunit e 1 Alternative name(s): D-serine-regulated transcript 1 protein Short name= DSR-1 V-ATPase 9.2 kDa membrane accessory protein V-ATPase M9.2 subu
Expression Region	2-81
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.