



# Recombinant Dog V-type proton ATPase subunit e 1 (ATP6V0E1)

<b>Product Code</b>	CSB-CF858984DO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9BDP4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Canis lupus familiaris (Dog) (Canis familiaris)
<b>Sequence</b>	AYHGLTVPLIVMSVFWGFGFCVPWFIPKGPNRGVIITMLVTCSVCCYLFWLIAI LAQLN PLFGPQLKNETIWYLKYHWP
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ATP6V0E1
<b>Protein Names</b>	Recommended name: V-type proton ATPase subunit e 1 Short name= V-ATPase subunit e 1 Alternative name(s): V-ATPase 9.2 kDa membrane accessory protein V-ATPase M9.2 subunit Vacuolar proton pump subunit e 1
<b>Expression Region</b>	2-81
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.