



Recombinant UPF0057 membrane protein ZK632.10 (ZK632.10)

Product Code	CSB-CF327352CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34655
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MCQILLAILAIFLPPIAVLLDVGCNCDLLINILLTCLGIIPGIIHAWYIILCKEKTVVQN IYVQTNDHGTAPPAYSPYSA
Source	in vitro E.coli expression system
Target Names	ZK632.10
Protein Names	Recommended name: UPF0057 membrane protein ZK632.10
Expression Region	1-80
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.