



Recombinant *Gossypium hirsutum* ATP synthase subunit c, chloroplastic (atpH)

Product Code	CSB-CF645685GHB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2L8Z2
Product Type	Transmembrane Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Sequence	MNPLISAASVIAAGLAVGLASIGPGVGQGTAAGQAVEGIARQPEAEGKIRGTL LSLAFM EALTIYGLVVALALLFANPFV
Source	in vitro <i>E.coli</i> expression system
Target Names	atpH
Protein Names	Recommended name: ATP synthase subunit c, chloroplastic Alternative name(s): ATP synthase F(0) sector subunit c ATPase subunit III F-type ATPase subunit c Short name= F-ATPase subunit c Lipid-binding protein
Expression Region	1-81
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.