



Recombinant Synechococcus elongatus ATP synthase subunit c (atpE)

Product Code	CSB-CF657433FPY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q31RF5
Product Type	Transmembrane Protein
Immunogen Species	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Sequence	MDSL TSAASVLAALAVGLAAIGPGIGQGSAAAGQAVEGIARQPEAEGKIRGTLL LSLAFM EALTIYGLVVALVLLFANPFA
Source	in vitro E.coli expression system
Target Names	atpE
Protein Names	Recommended name: ATP synthase subunit c Alternative name(s): ATP synthase F(0) sector subunit c F-type ATPase subunit c Short name= F-ATPase subunit c Lipid-binding protein
Expression Region	1-81
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.