



Recombinant Schizosaccharomyces pombe Succinate dehydrogenase cytochrome B subunit, mitochondrial (sdh3)

Product Code	CSB-CF525620SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O74882
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MFATRSFCLSSSLFRPAAQLLRPAGRSTLRNVWRRRSIATEHLTQTEANSRLAS QRVHRPN SPHLTIYEPQLTWYLSLHRIT
Source	in vitro E.coli expression system
Target Names	sdh3
Protein Names	Recommended name: Succinate dehydrogenase cytochrome B subunit, mitochondrial
Expression Region	1-82
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.