



Recombinant Rat Cytochrome b-c1 complex subunit 8 (Uqcrq)

Product Code	CSB-CF763178RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TQ16
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MGREFGNLTRIRHVISYLSLSPFEQRAFPHYFSKGI PNVLRRTRERILRVAPPFVLFYLIY TWGNQEFAQSKRKNPAKYENDK
Source	in vitro E.coli expression system
Target Names	Uqcrq
Protein Names	Recommended name: Cytochrome b-c1 complex subunit 8 Alternative name(s): Complex III subunit 8 Complex III subunit VIII Low molecular mass ubiquinone-binding protein Ubiquinol-cytochrome c reductase complex 9.5 kDa protein Ubiquino
Expression Region	1-82
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.