



Recombinant Pongo abelii NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 (NDUFA3)

Product Code	CSB-CF612636PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQ94
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	MAARVGAFRLRNTWDKEPVLVVSFVIGGLAVILPPLSPYFKYSIMINKATPYNYPV PVRDD GNMPDMPSHPQDPQGGSLEWLKKL
Source	in vitro E.coli expression system
Target Names	NDUFA3
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 Alternative name(s): Complex I-B9 Short name= CI-B9 NADH-ubiquinone oxidoreductase B9 subunit
Expression Region	1-84
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.