



Recombinant Human Amphiregulin (AREG)

Product Code	CSB-CF001986HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15514
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SVRVEQVVKPPQNKTESENTSDKPKRKKKGGKNGKNRRNRKKKNPCNAEFQ NFCIHGECK YIEHLEAVTCKCQQEYFGERCGEK
Source	in vitro E.coli expression system
Target Names	AREG
Protein Names	Recommended name: Amphiregulin Short name= AR Alternative name(s): Colorectum cell-derived growth factor Short name= CRDGF
Expression Region	101-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.