



Recombinant Human Protein TRIQK (C8orf83)

Product Code	CSB-CF730816HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q629K1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGRKDAATIKLPVDQYRKQIGKQDYKKTTPILRATKCLKAEAKKTAIGIKEVGLVLAAILA LLLAFYAFFYLRLTTDVPDLQDED
Source	in vitro E.coli expression system
Target Names	TRIQK
Protein Names	Recommended name: Protein TRIQK Alternative name(s): Triple repetitive-sequence of QXXK/R protein homolog
Expression Region	1-86
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.