



Recombinant Drosophila melanogaster Protein anon-73B1 (anon-73B1)

Product Code	CSB-CF524937DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O77286
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MSASADSLGAAAALDKYGDEDIFSLIRYGLYVGALFQFVCISAAVLMENNPDG QSNPES GEVTEREGEPVRTRLHKIRKLEKKRR
Source	in vitro E.coli expression system
Target Names	anon-73B1
Protein Names	Recommended name: Protein anon-73B1
Expression Region	1-87
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.