



Recombinant Uncharacterized 9.5 kDa protein in ORF3 5'region (xpaL2)

Product Code	CSB-CF859916BQT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q99165
Product Type	Transmembrane Protein
Immunogen Species	Bacillus licheniformis
Sequence	MKHIDKGTVVRTVLLFIALVNQTLIMFGKPVLPVAEDQIHTLADALYSAGSAAFTI AASL VAWYKNNYVTSKGKMQKEVLQKKGLTK
Source	in vitro E.coli expression system
Target Names	xpaL2
Protein Names	Recommended name: Uncharacterized 9.5 kDa protein in ORF3 5'region
Expression Region	1-87
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.