



Recombinant Uncharacterized protein ZK1098.9 (ZK1098.9)

Product Code	CSB-CF339475CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34608
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSEKHFILPSSMLMIVSAVFGGIGIITTIVFVILTVLHSKSAVCKPAGKEDMKKLN GIEG MQTIKEECGGSTETSSSKPKKKAKKEV
Source	in vitro E.coli expression system
Target Names	ZK1098.9
Protein Names	Recommended name: Uncharacterized protein ZK1098.9
Expression Region	1-87
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.