



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YCR024C-B (YCR024C-B)

Product Code	CSB-CF660350SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3E7Z8
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MCVCAIPFFEFFLPPFIPHYAFLLFVSSVRFTVNERCYLVCVLKLNCAFFFMVMI FELKR VCVSYLDRSRKIQIVSFFPFITIIFFHS
Source	in vitro E.coli expression system
Target Names	YCR024C-B
Protein Names	Recommended name: Uncharacterized protein YCR024C-B
Expression Region	1-88
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.