



Recombinant Salmonella typhimurium Secretion system apparatus protein SsaS (ssaS)

Product Code	CSB-CF302692SXB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P74891
Product Type	Transmembrane Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Sequence	MNDSELTQFVTQLLWIVLFTSMPVVLVASVVGIVSLVQALTQIQDQTLQFMIKL LAI AI TLMVSY PWLSGILLNYTRQIMLRIGEHG
Source	in vitro E.coli expression system
Target Names	ssaS
Protein Names	Recommended name: Secretion system apparatus protein SsaS
Expression Region	1-88
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.