



Recombinant Flagellar biosynthetic protein FliQ (fliQ)

Product Code	CSB-CF364731EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AC09
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MTPEsvMMMGTEAMKVALALAApLLLVALVTGLIISILQAATQINEMTLsFIPKIIA VFI AIIIAGPWMLNLLLLDYVRTLFTNLpYIIG
Source	in vitro E.coli expression system
Target Names	fliQ
Protein Names	Recommended name: Flagellar biosynthetic protein FliQ
Expression Region	1-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.