



Recombinant Bovine Short transient receptor potential channel 6 (TRPC6)

Product Code	CSB-CF865031BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9MYW0
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	IFIMVFVAFMIGMFNLYSYYIGAKQNEAFTTVEESFKTLFWAIFGLSEVKSVVINY NHKF IENIGYVLYGVYNVTMIVVLLNMLIAMIN
Source	in vitro E.coli expression system
Target Names	TRPC6
Protein Names	Recommended name: Short transient receptor potential channel 6 Short name= TrpC6
Expression Region	1-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.