



Recombinant Schizosaccharomyces pombe Dolichol-phosphate mannosyltransferase subunit 3 (dpm3)

Product Code	CSB-CF007136SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94633
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MQRHKKVILYYVSLTILYRVTYLFDLEEPWSTLRPYTPYLFILAFGSYLGITLLYNV ATT NDKPEAYVDLVKDIKEAQDALRSKGMTIED
Source	in vitro E.coli expression system
Target Names	dpm3
Protein Names	Recommended name: Dolichol-phosphate mannosyltransferase subunit 3 Alternative name(s): DPM synthase complex subunit 3 Dolichol-phosphate mannose synthase subunit 3 Dolichyl-phosphate beta-D-mannosyltransferase subunit 3 Mannose-P-doli
Expression Region	1-90
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.