



Recombinant Enterobacteria phage PRD1 Protein P18 (XVIII)

Product Code	CSB-CF335290EDO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27389
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage PRD1 (Bacteriophage PRD1)
Sequence	MPFGLIVIGIILAI AAYRDTLGELFSIIKDVSKDAKGFGYWVLA AVILGFAASIKPIK EP VNAFMILLMIVLLIRKRG AIDQISNQLRGS
Source	in vitro E.coli expression system
Target Names	XVIII
Protein Names	Recommended name: Protein P18 Alternative name(s): GpM Protein M
Expression Region	1-90
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.