



Recombinant Rat Activator of apoptosis harakiri (Hrk)

Product Code	CSB-CF010741RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P62817
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MCPCPRHRGRGPPAVCGCGDARPGLRWAAAQVTALRLQALGDELHRRAMR RRARPRDPLP ALLPALRARWPWLCAAQAALAAWLLGRRSA
Source	in vitro E.coli expression system
Target Names	Hrk
Protein Names	Recommended name: Activator of apoptosis harakiri Alternative name(s): BH3-interacting domain-containing protein 3 Neuronal death protein DP5
Expression Region	1-92
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.